



## PRODUCT DATA SHEET

### HISTOWAX POLYMER ADDED WITH DMSO (PHC9010)

A highly purified paraffin wax with added synthetic polymers and a small addition of DMSO for faster tissue infiltration. Suitable for all general embedding and sectioning use. Characterised by low shrinkage on cooling, excellent sectioning performance with no expansion. Superb cutting ease and ribbon continuity. Ideal for larger tissue samples.

<b>Congealing point</b>	<b>54 - 57°C</b>
<b>Melt Range</b>	<b>56 – 58°C</b>
<b>Visual</b>	<b>White free flowing pellets</b>
<b>Flocculation (48hrs @ 63.0°C)</b>	<b>Complies with PHC flocculation test.</b>

Available in 1Kg, 2Kg, 2.5Kg bags as well as 10Kg and 25Kg cartons.

*This is a general guide for reference purposes only. Please contact us quoting the PHC reference number for detailed specifications.*

